新冠抗原检测试剂产品手册
**Product Manual of
COVID-19 Antigen Detection Kit**

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COMPANY PROFILE

New Gene (Hangzhou) Bioengineering Co., Ltd. is located in Hangzhou, China. It is a high-tech company engaged in the research, development, manufacture and distribution of biological products. It is committed to creating biological materials such as antigens and antibodies, in vitro diagnostic reagents and related devices, and also the complete industrial chain of artificial intelligence assisted diagnosis system. The product line covers a full range of in vitro diagnostic products such as immune diagnosis, molecular diagnosis, and microbiological testing. NEWGENE has profound technical accumulation and unique technological advantages in the areas of early cancer screening, rapid detection of infectious diseases, and rapid screening of geriatric diseases.

NEWGENE's manufacturing system meets GMP standards for medical devices, and is certified with ISO13485 by British BSI. Relevant in vitro diagnostic reagent products have obtained the EU CE certification. NEWGENE is also a member on the "allow list" issued by Chinese Ministry of Commerce for anti-epidemic products exporting.

At present, NEWGENE COVID-19 Antigen Detection Kit has **registered in** many countries, including **Germany, France, Italy, Switzerland, Belgium, Portugal, Czech, Denmark, Hungary, Greece, Poland, Sweden, Moldova, Peru, Argentina, Ecuador, Kenya, Zimbabwe, South Africa, Malaysia, Thailand, Brunei, Congo etc.**, and passed the clinical verification of multiple national laboratories. The products are suitable for children under 14 years old.

Self-test Approval of EU

CERTIFICATE

EC Certificate No. 1434-IVDD-449/2021**EC Design-examination
Directive 98/79/EC concerning
in vitro diagnostic medical devices**Polish Centre for Testing and Certification certifies
that manufactured by:**New Gene (Hangzhou) Bioengineering Co., Ltd.
Room 1606, Floor 16, Building 5, 688 Bin'an Road, Changhe Street,
Binjiang District, Hangzhou City, Zhejiang Province,
P. R. China**in vitro diagnostic medical devices
for self-testing**COVID-19 Antigen Detection Kit - Nasal Swab**in terms of design documentation, comply with requirements
of Annex III (Section 6) to Directive 98/79/EC (as amended)
implemented into Polish law,

as evidenced by the audit conducted by the PCBC

Validity of the Certificate: from 11.08.2021 to 27.05.2024

The date of issue of the Certificate: 11.08.2021

The date of the first issue of the Certificate: 11.08.2021

Issued under the Contract No. MD-116
Application No: 239/2021
Certificate bears the qualified signature.
Warsaw, 11.08.2021
Module A1
Elektronicznie
podpisany przez Anna
Małgorzata Wyroba
Data: 2021.08.11
09:14:18 +02'00'
Vice-President

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Manufacturer: New Gene (Hangzhou) Bioengineering Co., Ltd.
Address: Room 1606, Floor 16, Building 5, 688 Bin'an Road, Changhe Street,
Binjiang District, Hangzhou City, Zhejiang Province, P. R. China

EC Representative: SUNGO Europe B.V.
Address: Olympisch Stadion 24, 1076DE Amsterdam, Netherlands
EC Certificate No.: 1434-IVDD-449/2021

Product Name: COVID-19 Antigen Detection Kit – Nasal Swab
Specification: 1Test/Box, 5Tests/Box, 25Tests/Box
Classification: Self Test (IVDD)
Conformity Assessment
Procedure: Annex III (Section 6) to Directive 98/79/EC

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Signature: 
Name/ Position: Mingfu Li / General Manager
Date: 11/08/2021
Place: Hangzhou, Zhejiang, China

By Royal Charter

Certificate of Registration

QUALITY MANAGEMENT SYSTEM - ISO 13485:2016

This is to certify that: New Gene (Hangzhou)
Bioengineering Co., Ltd.
Room 1606, 16th Floor, No.5 Building
688 Bin'an Road
Binjiang District
Hangzhou
Zhejiang
310052
China

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Holds Certificate No: **MD 729179**

and operates a Quality Management System which complies with the requirements of ISO 13485:2016 for the following scope:

Design and Development, Manufacture and Distribution of In-vitro Diagnostic Rapid Test Kit of Drug Abuse, Manufacture and Distribution of In-vitro Diagnostic Rapid Test Kit of Infectious Diseases.

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For and on behalf of BSI:

Gary E Slack, Senior Vice President - Medical Devices

Original Registration Date: 2020-07-27

Latest Revision Date: 2020-07-27

Effective Date: 2020-07-27

Expiry Date: 2023-07-26

Page: 1 of 1

...making excellence a habit.™

Self-Testing

CE 1434

COVID-19 Antigen Detection Kit - Nasal Swab

No.	Components	25 Tests/Box	5 Tests/Box	1 Test/Box
1	Test Card	25	5	1
2	Sample Extraction Tube & Tube Cap	25	5	1
3	Sampling Swab: <i>for Nasal Swab</i>	25	5	1
4	Package Insert	1	5	1

25 Tests/Box


5 Tests/Box


1 Test/Box

COMPONENTS


Aluminum Foil Pouch


Test Card
(Packed in the pouch)


Desiccant
(Discard it.
Do not open)

**Sample Extraction Tube &
Tube Cap**

Sampling Swab

Instructions for Use
COVID-19 Antigen Detection Kit
EN

For self-test use / Suitable for non-professionals to conduct self-test.

PRECAUTIONS BEFORE USING THE PRODUCT


1. Read the instructions carefully prior to first use.
2. For people who are not able to perform the test themselves, the test should be conducted by the legal guardians.
3. For children under the age of 15, the self-test should be conducted under adult supervision.
4. This test detects SARS-CoV-2 antigen in nasal cavity secretions, which is collected by a sterile nasal swab.
5. For people who has recent nasal trauma or surgery, or has severe coagulopathy, gentle operation is required for nasal swab collection to avoid injuries to the nose.
6. Please use the components provided in the kit for testing. Do not use components from other sources.
7. Please use this product in a place with sufficient light, so as to interpret the results accurately.

PRECAUTIONS AFTER USING THE PRODUCT


1. If you get a positive result, please contact your family doctor, or seek help from a professional medical facility as soon as possible. You need a nucleic acid test to confirm viral infection.
2. A negative result cannot completely exclude the possibility of viral infection. Incorrect sampling or low viral load may also cause a false negative result.

Instructions for Use

TEST PROCEDURES


INTERPRETATION OF RESULTS


Positive (+): Red bands appear at both of T and C line in 15 to 30 minutes.

A white band at the T line should be considered as a negative result.

Negative (-): A red band appears at C line while no red band appears at T line in 15 to 30 minutes after sample loading.

Invalid: If no red band appears at C line, it indicates that the test result is invalid.

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
Dr. M.J. van de Velde

DECLARATION OF CONFORMITY

Regarding In Vitro Diagnostic Directive (98/79/EC)

Manufacturer: New Gene (Hangzhou) Bioengineering Co., Ltd.
Address: Room 1606, Floor 16, Building 5, 688 Bin'an Road, Changhe Street, Binjiang District, Hangzhou City, Zhejiang Province, P. R. China

EC Representative: SUNGO Europe B.V.
Address: Olympisch Stadion 24, 1076DE Amsterdam, Netherlands

Product Name: COVID-19 Antigen Detection Kit
Product Code: COVID-19-NG08
Specification: 25Tests/Box 1Test/Box
Classification: Others (IVDD)

Conformity Assessment Procedure: Annex III of In Vitro Diagnostic Directive (98/79/EC)

We herewith declare that the above-mentioned products meet the requirements of In Vitro Diagnostic Directive (98/79/EC) and the following harmonized standards.

- | | |
|-------------------|---------------------|
| EN 23640-2015 | EN 13640:2002 |
| EN 980:2016 | EN 13641:2002 |
| EN ISO 14971:2019 | EN ISO 18113-1:2011 |
| EN 13612:2002 | EN ISO 18113-4:2011 |

On behalf of SUNGO Europe office, I confirmed we are EU REP of the company who issue this document.

Signature: 
Name/ Position: Mingfu Li / General Manager
Date: 29/09/2020
Place: Hangzhou, Zhejiang, China

 Authorized Signature (S)

COVID-19 Antigen Detection Kit (Nasal Swab)

For Professional Use

No.	Components	25 Tests/Box	5 Tests/Box	1 Test/Box
1	Test Card	25	5	1
2	Sample Extraction Tube & Tube Cap	25	5	1
3	Sampling Swab: for Nasal Swab	25	5	1
4	Package Insert	1	5	1

25 Tests/Box


5 Tests/Box


1 Test/Box


Test Procedures: Sampling

Test Procedures: Detection

Interpretation of Results

Positive(+): Red bands appear at both of T and C line in 15 to 30 minutes. A white band at the T line should be considered as a negative result.

Negative(-): A red band appears at C line while no red band appears at T line in 15 to 30 minutes after sample loading.

Invalid: As long as no red band appears at C line, it indicates that the test result is invalid, and should retest the sample with another test card.

Performance

Sensitivity	Specificity
97.1%	99.2%

COVID-19 Antigen Detection Kit (Nasopharyngeal Swab)

For Professional Use

No.	Components	25 Tests/Box	5 Tests/Box	1 Test/Box
1	Test Card	25	5	1
2	Sample Extraction Tube & Tube Cap	25	5	1
3	Sampling Swab: <i>for Nasopharyngeal Swab</i>	25	5	1
4	Package Insert	1	5	1

25 Tests/Box


5 Tests/Box


1 Test/Box


Test Procedures: Sampling

Test Procedures: Detection

Interpretation of Results

Positive(+): Red bands appear at both of T and C line in 15 to 30 minutes. A white band at the T line should be considered as a negative result.

Negative(-): A red band appears at C line while no red band appears at T line in 15 to 30 minutes after sample loading.

Invalid: As long as no red band appears at C line, it indicates that the test result is invalid, and should retest the sample with another test card.

Performance

Sensitivity	Specificity
98.0%	99.1%

COVID-19 Antigen Detection Kit (Oropharyngeal Swab)

For Professional Use

No.	Components	25 Tests/Box	5 Tests/Box	1 Test/Box
1	Test Card	25	5	1
2	Sample Extraction Tube & Tube Cap	25	5	1
3	Sampling Swab: for Oropharyngeal Swab	25	5	1
4	Package Insert	1	5	1

25 Tests/Box


5 Tests/Box


1 Test/Box


Test Procedures: Sampling

Test Procedures: Detection

Interpretation of Results

Positive(+): Red bands appear at both of T and C line in 15 to 30 minutes. A white band at the T line should be considered as a negative result.

Negative(-): A red band appears at C line while no red band appears at T line in 15 to 30 minutes after sample loading.

Invalid: As long as no red band appears at C line, it indicates that the test result is invalid, and should retest the sample with another test card.

Performance

Sensitivity	Specificity
95.7%	99.0%

COVID-19 Antigen Detection Kit (Saliva)

For Professional Use

No.	Components	25 Tests/Box	5 Tests/Box	1 Test/Box
1	Test Card	25	5	1
2	Sample Extraction Tube & Tube Cap	25	5	1
3	Paper Cup	25	5	1
4	Sputum Dropper	25	5	1
5	Package Insert	1	5	1

25 Tests/Box


5 Tests/Box


1 Test/Box


Test Procedures: Sampling

Test Procedures: Detection

Interpretation of Results

Positive(+): Red bands appear at both of T and C line in 15 to 30 minutes. A white band at the T line should be considered as a negative result.

Negative(-): A red band appears at C line while no red band appears at T line in 15 to 30 minutes after sample loading.

Invalid: As long as no red band appears at C line, it indicates that the test result is invalid, and should retest the sample with another test card.

Performance

Sensitivity	Specificity
97.3%	99.2%

PACKAGING INFORMATION For Professional Use

Nasopharyngeal Swab: **NPS**

Nasal Swab: **NS**

Oropharyngeal Swab: **OS**

Saliva: **S**

25 Tests/Box

Sample	NPS	NPS+S	NS	NS+S	OS	OS+S	S
Box (mm)	230*140*80						230*120*67
Box weight (kg)	0.36						
Carton (mm)	585*485*425						510*490*360
Carton weight (kg)	1.5						1.3
PCS/Box	25						
Boxes/Carton	40						
PCS/Carton	1000						
Volume/Carton	0.12CBM						0.09CBM
NW/Carton (kg)	14.3						15.6
GW/Carton (kg)	15.8						16.9

5 Tests/Box

Sample	S	NS	NS+S	NPS	OS	NPS+S	OS+S
Inner box (mm)	193*85*42						
Outer box (mm)	225*197*89						
Outer box weight (kg)	0.5						
Carton (mm)	470*410*470						
Carton weight (kg)	1.4						
PCS/Inner Boxes	5						
Inner Boxes/Outer Box	20						
PCS/Carton	500						
Volume/Carton	0.09CBM						
NW/Carton (kg)	10						
GW/Carton (kg)	11.4						

1 Test/Box

Sample	S	NS	NS+S	NPS	OS	NPS+S	OS+S
Inner box (mm)	143*83*15			170*66*15			
Inner box weight (kg)	0.027			0.028			
Outer box (mm)	305*197*88			277*182*112			
Outer box weight (kg)	0.8			0.83			
Carton (mm)	630*420*470			590*570*395			
Carton weight (kg)	1.8			2.2			
PCS/Inner Boxes	25						
Inner Boxes/Outer Box	20						
PCS/Carton	500						
Volume/Carton	0.13CBM			0.133CBM			
NW/Carton (kg)	16.5			16.6			
GW/Carton (kg)	18.3			18.8			

PACKAGING INFORMATION Self-Testing

25 Tests / Box

Sample	Nasal Swab
Box (mm)	230x120x75
Box weight (kg)	0.27
Carton (mm)	520X480X400
Carton weight (kg)	1.5
PCS/Box	25
Boxes/Carton	40
PCS/Carton	1000
Volume/Carton	0.1CBM
NW/Carton (kg)	11
GW/Carton (kg)	12.5

5 Tests / Box

Sample	Nasal Swab
Inner box (mm)	122x68x49
Outer box (mm)	258X128X75
Outer box weight (kg)	0.38
Carton (mm)	550X530X395
Carton weight (kg)	1.4
PCS/Inner Boxes	5
Inner Boxes/Outer Box	40
PCS/Carton	1000
Volume/Carton	0.11CBM
NW/Carton (kg)	15.6
GW/Carton (kg)	17kg

1 Test / Box

Sample	Nasal Swab
Inner box (mm)	122X83X15
Outer box (mm)	255X170X105
Outer box weight (kg)	0.68
Carton (mm)	545X530X370
Carton weight (kg)	1.8
PCS/Inner Boxes	25
Inner Boxes/Outer Box	20
PCS/Carton	500
Volume/Carton	0.10CBM
NW/Carton (kg)	12.45
GW/Carton (kg)	14.25

EUROPEAN COMMISSION
DIRECTORATE-GENERAL FOR HEALTH AND FOOD SAFETY

Public health, country knowledge, crisis management
Health Security

EU health preparedness: A common list of COVID-19 rapid antigen tests and a common standardised set of data to be included in COVID-19 test result certificates

Agreed by the Health Security Committee

This document was agreed by the HSC on 17 February 2021

Annex I

Common list of COVID-19 rapid antigen tests

A first update was agreed by the HSC on 10 May 2021; A second update was agreed by the HSC on 16 June 2021; A third update was agreed by the HSC on 7 July 2021; A fourth update was agreed by the HSC on 14 July 2021; A fifth update was agreed by the HSC on 23 July 2021.

IMPORTANT: A (interim) grace period of 8 weeks applies whenever updates are made to Annex I, the common list of COVID-19 rapid antigen tests

New Gene (Hangzhou) Bioengineering Co., Ltd.	COVID-19 Antigen Detection Kit	Yes	98% sensitivity Nasal swab	DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 92,5% at <Ct30 and 100% at <Ct25)		DE ^[2]		DE ^[2]		1501	16 June 2021
----------------------------------------------	--------------------------------	-----	-------------------------------	-------------------------------------------------------------------------------------------------------	--	-------------------	--	-------------------	--	------	--------------

Validation protocol for the Newgene bioengineering COVID-19 Rapid Antigen Test

Testing laboratory:
Molecular Diagnostic Laboratory
Via Petrini 2
CH-6900 Lugano
Switzerland

Discussion and Conclusion.

In conclusion, the Newgene antigen detection test is highly precise and accurate (**100% specificity and 95.1% sensitivity**), it is non-invasive and mimics the PCR results very closely where PCR is considered the golden standard technique. The device can easily be operated and used by non-medically trained personnel, does not need a laboratory setting and could be intended for regular use by regular people.

Lugano, June 18th 2021.
Dr. G. Soldati
CEO
Molecular Diagnostic Laboratory
Via Petrini 2
CH-6900 Lugano, Switzerland

Evaluation Report (Children under the age of 14)

Discussion

In the present study, the COVID-19 Antigen Detection Kit - Nasal Swab by New Gene (Hangzhou) Bioengineering Co., Ltd. has shown highly reliable performance in sample from children under the age of 14. Compared to conventional RT-PCR tests, the rapid antigen tests can meet a wider range of test needs. As the antigen tests take only 20 to 30 minutes, they are feasible for use in emergency scenarios where a test result is demanded immediately. Also, the rapid antigen tests doesn't require special instruments and training to use, they are also capable for resource limited scenarios like point-of-care testing and self-test by laypeople. Therefore, the implementation of rapid antigen test may totally change the strategies to control COVID-19. Community residents can conduct the rapid antigen test in a frequent manner, like twice or three times a week, to identify COVID-19 cases in the early stage of infection. This strategy may help to stop the transmission of COVID-19 as early as possible. In summary, the COVID-19 Antigen Detection Kit - Nasal Swab has shown satisfying sensitivity, specificity, and total accuracy in the present evaluation. It can be used as a rapid tool to assist the early diagnosis of COVID-19 cases in children under the age of 14.

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Evaluation Report (Saliva)

Result Analysis:

The test results of RT-PCR and the COVID-19 Antigen Detection Kit are summarized in 2x2 table below.

		RT-PCR Reagent		Total
		Positive	Negative	
Antigen Test	Positive	204	0	204
	Negative	16	200	216
Total		220	200	420

Sensitivity (%) = $204 \div (204 + 16) \times 100\% = 92,73\%$; 95% CI: 88,46% ~ 95,79%

Specificity (%) = $200 \div (0 + 200) \times 100\% = 100,0\%$; 95% CI: 98,17% ~ 100,00%

Total accuracy (%) = $(204 + 200) \div (204 + 0 + 16 + 200) \times 100\% = 96,19\%$;

95% CI: 93,89% ~ 97,81%

The 95% Confidence Intervals of sensitivity, specificity, and total accuracy are calculated following the binomial distribution.

Analysis results by Kappa consistency test shows that the Kappa value = 0,9239 (95% CI: 0,8875 ~ 0,9604). As $Kappa \geq 0.75$, it suggests good consistency between the COVID-19 Antigen Detection Kit and the COVID-19 RT-PCR detection reagent. Therefore, the COVID-19 Antigen Detection Kit has shown reliable performance on detecting SARS-CoV-2 virus in saliva samples from suspected COVID-19 patients.

(+++) - positive - intensity of line T similar or higher than line C

(++) - positive - intensity of line T lower than line C

(+) - positive - border intensity of line T

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
Registration or Allowed List (Partially)

Self-test Approval in Malaysia

Home / Announcement / SELF-TEST COVID-19 TEST KIT FOR CONDITIONAL APPROVAL (APPROVED)

SELF-TEST COVID-19 TEST KIT FOR CONDITIONAL APPROVAL (APPROVED)

SELF-TEST COVID-19 TEST KIT FOR CONDITIONAL APPROVAL (APPROVED)

The list of Self-Test Covid-19 Test Kit that is approved for Conditional Approval based on the decision on the consensus of the Covid-19 Test Kit Expert committee is as follows.

All test submissions are scored according to: -

- the manufacturer reported clinical and analytical performance evidence.
- the evaluation test results from testing facilities are according to the committee evaluation criteria set by Clinical expert panels.
- Supporting Documents for COVID-19 IVD Test Kits Conditional Approval.

The use of COVID-19 self- test kit shall be limited for screening purpose only and all test results need further confirmation using RT-PCR.

This test kit can be supplied by registered pharmacists or private healthcare facilities.

Below is the list of all tests that have been selected to date and the status is Conditional Approval (please note: list is updated on a routine basis).

Notes:

*Sample type is based on testing facility evaluation report.

NO	COMPANY NAME	PRODUCT NAME	MANUFACTURER	IDENTIFIER	DETECTION	SAMPLE TYPE
11	Dewina Consult Sdn Bhd	NEWGENE Bioengineering COVID-19 Antigen Detection Kit	New Gene (Hangzhou) Bioengineering Co., Ltd. P.R. China	COVID-19-NG08	RTK-Antigen (Self-test)	Saliva or Nasal swab

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A gyártó neve: New Gene (Hangzhou) Bioengineering Co.Ltd.
A meghatalmazott képviselő neve: Sungo Europe B.V.